# **RIDASCREEN<sup>®</sup> Saxitoxin**

Enzymimmunoassay zur quantitativen Bestimmung von  
Saxitoxin (Paralytic Shellfish Poison, PSP)

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# RIDASCREEN® Saxitoxin

Enzyme immunoassay for the quantitative analysis of saxitoxin (Paralytic Shellfish Poison, PSP)

## Brief information

RIDASCREEN® Saxitoxin (Art. No.: R1901)

Competitive enzyme immunoassay for the quantitative analysis of saxitoxin and related toxins in mussels.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: washing, homogenization, boiling with HCl and centrifugation

Time requirement: sample preparation (for 10 samples) ....approx. 1.5 h  
test implementation .....1.5 h  
(regardless of the number of samples)

Detection limit: 2 ppb

Recovery rate: approx. 90 %

Cross-reactivity: Saxitoxin ..... 100 %  
Decarbamoyl saxitoxin ..... 10 - 30 %  
Gonyautoxins II, III, B<sub>1</sub>, C<sub>1</sub> und C<sub>2</sub>..... 10 - 30 %

The specificity of the RIDASCREEN® Saxitoxin assay was established by analyzing the cross reactivity to the corresponding toxins. Because of the fact, that pure standard substances are not available, a contamination of the corresponding toxins with saxitoxin is possible. Therefore, the results obtained indicate approx. the above mentioned or less cross reactivities.

## 1. Intended use

The RIDASCREEN® Saxitoxin test is a competitive enzyme immunoassay for the quantitative analysis of saxitoxin and related toxins in mussels.

## 2. General

Saxitoxin is a potential neurotoxin which occurs in various algae (Gonyaulax-species). Over the food chain it could reach mussels. Therefore, saxitoxin is then mainly responsible for the high toxicity of some mussels.

Saxitoxin, the principal component of paralytic shellfish poisons (PSP), is responsible for respiratory paralysis which in 8 % of cases leads to death. Saxitoxin blocks the sodium ion channels in nerve and muscle membranes. The lethal dose in humans is 1 to 3 mg. Numbness and respiratory arrest can occur already after oral ingestion of as little as 0.5 – 1.0 µg saxitoxin.

An established and reliable detection method for these algae toxins has been the "mouse mean death time bioassay". The limitations of this procedure are the high variability of the results and low sensitivity (approx. 37 µg PSP/100 g of sample). This low sensitivity is particularly problematic given the allowable limits of 40 and 80 µg PSP/100 g of sample respectively set by the individual member states of the EU. In addition, this method is to be critically reconsidered from the animal protection point of view.

As an alternative, the RIDASCREEN® Saxitoxin assay was developed for the rapid and sensitive detection of Saxitoxin and related toxins in mussels.

*In a comparative study using over 30 different naturally contaminated mussel samples the RIDASCREEN® Saxitoxin assay showed excellent correlation with the fluorimetric method (the german official method according to §35 LMBG).*

## 3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with antibodies directed against saxitoxin. Standards or sample solutions and enzyme conjugated saxitoxin are added. Free and enzyme conjugated saxitoxin compete for the antibody binding sites (competitive enzyme immunoassay). Any unbound enzyme conjugate is then removed in a washing step. Enzyme substrate (urea peroxide) and chromogen (tetramethylbenzidine) are added to the wells and incubated. Bound enzyme conjugate converts the colorless chromogen into a blue product.

The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm (optional reference wavelength  $\geq 600$  nm). The absorption is inversely proportional to the saxitoxin concentration in the sample.

#### 4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)  
coated with antibodies directed against saxitoxin
- 6 x Standard solutions, 1.3 ml each  
0 ppt (zero standard), 10 ppt, 30 ppt, 90 ppt, 270 ppt, 810 ppt  
saxitoxin in aqueous solution, ready to use
- 1 x Conjugate (0.7 ml) .....red cap  
peroxidase conjugated saxitoxin, concentrate
- 1 x Substrate (7 ml solution).....green cap  
contains urea peroxide
- 1 x Chromogen (7 ml) ..... blue cap  
contains tetramethylbenzidine
- 1 x Stop solution (14 ml) .....yellow cap  
contains 1 N sulfuric acid
- 1 x Buffer (50 ml)  
Sample and conjugate dilution buffer

#### 5. Materials required but not provided

##### 5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- mixer
- magnetic stirrer
- centrifuge
- graduated pipettes
- 10  $\mu$ l-, 50  $\mu$ l-, 100  $\mu$ l- and 500  $\mu$ l-micropipettes

## 5.2. Reagents:

– 0.1 M HCl (hydrochloric acid)

## 6. Warnings and precautions for the users

Standards contain saxitoxin. Particular care should be taken. Do not allow the reagent to come in contact with the skin (use gloves).

The stop solution contains 1 N sulfuric acid. Avoid contact of the reagent with the skin.

Do not use RIDASCREEN® Saxitoxin kit past the kit expiration date on the label. Dilution or adulteration of these reagents may result in loss of sensitivity.

Do not interchange individual reagents between kits of different lot numbers.

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). DO NOT FREEZE.

Return any unused microwells to their original foil bag and reseal them together with the desiccant provided.

The colorless chromogen is light sensitive, therefore, avoid exposure to direct light.

## 8. Indication of instability or deterioration of reagents

Any coloration of the chromogen solution is indicative of deterioration and the reagent should be discarded.

A value of less than 0.6 absorbance units ( $A_{450 \text{ nm}} < 0.6$ ) for the zero standard may indicate deterioration of reagents.

## 9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

### 9.1. mussels

- after removing the shell, wash the mussel meat with water and homogenize
- mix 10 g of homogenized mussels with 10 ml HCl (0.1 M) and boil for 5 min while stirring
- centrifuge: 10 min /approx. 3500 g / room temperature (20 - 25 °C / 68 - 77 °F) (possible clouding of the supernatant after centrifugation does not have an effect on the test result)
- use 20 µl of the supernatant and fill up to 2 ml with sample dilution buffer
- employ 50 µl per well in the assay

## 10. Test implementation

### 10.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. Return all reagents to 2 - 8 °C (36 - 46 °F) immediately after use.
3. Do not allow microwells to dry between working steps.
4. Reproducibility in any EIA is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.
5. Avoid direct sunlight during all incubations. Covering the microtiter plates is recommended.

### 10.2. Saxitoxin enzyme conjugate

The saxitoxin enzyme conjugate (bottle with red cap) is provided as a concentrate. Since the diluted enzyme conjugate has a limited stability, only the amount which actually is needed should be reconstituted. Before pipetting, the enzyme conjugate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 in buffer (e. g. 200 µl conjugate concentrate + 2 ml buffer, ready to use sufficient for 4 microtiter strips).

### 10.3. Antibody-coated microtiter strips

The foil bag is cut open along the transverse side beyond the zip. The required wells are extracted together with the frame. Those wells not required are kept together with the drying agent, well sealed in the foil bag and should continuously be stored at 2 - 8 °C (36 - 46 °F).

### 10.4. Test procedure

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50µl of standard or prepared sample to separate duplicate wells.
3. Add 50 µl of diluted enzyme conjugate solution to the bottom of each well. Mix gently by rocking the plate manually and incubate for 1 h at room temperature (20 - 25 °C / 68 - 77 °F).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all wells with 250 µl of distilled water and pour out the liquid again. Repeat two more times.
5. Add 50 µl of substrate and 50 µl of chromogen to each well. Mix gently by rocking the plate manually and incubate for 30 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µl of the stop solution to each well. Mix gently by rocking the plate manually and measure the absorbance at 450 nm against an air blank. Read within 60 minutes after addition of stop solution.

## 11. Results

The mean values of the absorbance values obtained for the standards and the samples are divided by the absorbance value of the first standard (zero standard) and multiplied by 100. The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages.

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the saxitoxin concentration in [ng/kg]. The calibration curve should be virtually linear in 30 - 270 ng/kg (ppt) range. The saxitoxin concentration in ng/kg corresponding to the absorbance of each sample can be read from the calibration curve.

In order to obtain the saxitoxin concentration in ng/kg actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

mussels ..... 200

**Please note:**

For evaluation of the ELISA kits special software for the RIDASCREEN® tests has been developed by R-Biopharm. The RIDA® SOFT Win can be ordered by your local distributor.

**12. Sensitivity**

The mean lower detection limit of the RIDASCREEN® Saxitoxin test is 10 ng/kg. According to the sample preparation protocol, the detection limit is approx. 2 µg/kg (ppb) saxitoxin in mussels.

Fig. 1: Calibration curve of a RIDASCREEN® Saxitoxin kit

### 13. Specificity

The specificity of the RIDASCREEN® Saxitoxin test was determined by analyzing the cross-reactivities to corresponding toxins.

Cross-reactivity:	Saxitoxin .....	100 %
	Decarbamoyl saxitoxin .....	10 - 30 %
	Gonyautoxins II, III, B <sub>1</sub> , C <sub>1</sub> und C <sub>2</sub> .....	10 - 30 %

The specificity of the RIDASCREEN® Saxitoxin assay was established by analyzing the cross reactivity to the corresponding toxins. Because of the fact, that pure standard substances are not available, a contamination of the corresponding toxins with saxitoxin is possible. Therefore, the results obtained indicate approx. the above mentioned or less cross reactivities.

### 14. Reproducibility

The precision within a series was determined from the results of three different experiments. Fig. 2 shows the interassay precision profile of RIDASCREEN® Saxitoxin. The coefficients of variation (% CV) obtained for the absorption values of the standards are entered against the corresponding saxitoxin concentrations. The coefficients of variation are so low with respect to the whole range of the figure that a high reproducibility of the results is ensured.


Abb. 2: Interassay precision profile of RIDASCREEN® Saxitoxin

## 15. Recovery rate

The recovery rate for artificial contaminated mussel samples is approx. 90 % with a mean coefficient of variation of 12 %.

## Appendix

### Literature

1. Usleber, E. et al.; Letters in Applied Microbiology, 13, 275-277 (1991):  
Direct enzyme immunoassay in microtitration plate and test strip format for the detection of saxitoxin in shellfish
2. Schantz, E. J.; American Chemical Society, 9, 99-111 (1984):  
Historical Perspective on Paralytic Shellfish Poison
3. Sullivan, J. J. and Iwaoka, W. T.; J. Assoc. Off. Anal. Chem., Vol. 66, 2, 297-303 (1983):  
Seafood Toxins: High Pressure Liquid Chromatographic Determination of Toxins Associated with Paralytic Shellfish Poisoning

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG